

Atraumatische Gefäßklemmen
 Atraumatic Vascular Clamps
 Pinces vasculaires atraumatiques
 Pinzas vasculares atraumáticas
 Pinze vascolari atraumatiche

0013-000F


ATRAUMA




DE BAKEY

Tangentialklemmen
 Tangential Occlusion Clamps
 Pinces-clamps pour occlusion tangentielle
 Pinzas para oclusión tangencial
 Pinze per occlusione tangenziale

*Jaw tips should close
 3th ratcheted
 completely close
 7-8th ratcheted*




SATINSKY



BD-523-25
 25 cm/10"

BD-523-24
 24 cm/9 1/2"

BD-523-20
 20 cm/8"



BD-523-26
 26 cm/10 1/4"

BD-523-27
 27 cm/10 1/2"

CM-Instrumente GmbH
 Gärtnacker 56
 D-78532 Tuttlingen / Germany
 Tel.: +49 (0) 74 62-204 20-70
 Fax: +49 (0) 74 62-204 20-71